PAL[®] Infiltration System

C



МАХ

MICROAIRE" PAPOSCULPTOR"

MIN

MICROAIRE

• PAL®

3

PAL[®] Infiltration System

The PAL[®] Infiltration Pump integrates with the PAL[®] System to activate the peristaltic pump from the PAL-650 Handpiece. Custom electronic controls ensure consistent flow at low speeds and a cooling fan prevents overheating of the motor from heavy use. PAL[®] Infiltration Tubing eliminates obstructions with a large-bore design and smooth connections.



System Accessories



Infiltration Pump to 5020 Connector Cable INF-CBL-5020



Infiltration Foot Switch INF-FOOT-1



PAL® Infiltration Tubing PAL-INF-1600



PAL® Single-Use Infiltration Cannulas PAL® Multi-Use Infiltration Cannulas PAL-INF-XXXXX PAL-INF-RXXXXX

A versatile system that expands the capabilities of PAL[®]



Ordering Information

PAL [®] Infiltration System		Qty.
INF-PUMP	PAL [®] Infiltration Pump	EA
INF-CBL-5020	Infiltration Pump to 5020 Connector Cable	EA
INF-FOOT-1	Infiltration Foot Switch	EA
PAL-INF-1600	PAL® Infiltration Tubing	5PK
PAL-INF-XXXXX	PAL [®] Single-Use Infiltration Cannulas	EA
PAL-INF-RXXXXX	PAL [®] Multi-Use Infiltration Cannulas	EA

To order, call toll free 800.722.0822

Learn more about our products, including the PAL® Infiltration System, at www.microaire.com.



MicroAire Surgical Instruments 3590 Grand Forks Boulevard Charlottesville, VA 22911 USA 800 722 0822 or 434 975 8000 www.microaire.com

© 2022 MicroAire Surgical Instruments, LLC All rights reserved. LIT-INF-PUMP Rev B

